

Tools and Test Equipment

1. KoKo System Diagnostics program.
2. Conditioning Syringe (K700953).
3. Silicone lubrication (K022772).
4. Phillips screwdrivers
5. Flat bladed screwdrivers
6. Dental style pick.
7. Allen Keys (Imperial)
8. RS Trimmer Tool
9. Contact cleaner (314957)
10. Dwyer 0 - 5 ltr Flow Meter
11. Comark Pressure Meter (Or equivalent)
12. Calibrated Manometer (10 cmH2O)



Specialist Tools You Will Need